Optical Clearing Innovation

www.sunjinlab.com **Product Information**

iSpacer® (imaging Spacer)

An imaging spacer keeping the biological tissue in good shape

INTRODUCTION

iSpacers are made from adhesive tapes in different thickness. Simply press a spacer to a microscope slide or coverslip, a sealed watertight well is formed to contain RapiClear in place and prevents evaporation. The well is sturdy and provide supreme support to thick and free-floating specimens, allowing preservation of tissues' delicate internal structures without compression. The specimen and spacer can be sandwiched between two coverslips. To better fulfil the user's different application needs, iSpacers not only come in various thicknesses, there are also iSpacers with single-sided stickers, as well as double-sided stickers. Utilizing these various designs, iSpacers can be stacked to create different depths desired for direct use with confocal microscopy.

OPERATION


Single-sided sticky iSpacers

- 1. Clean the microscope slide and coverslip (size: 22 x 22 mm, 24 x 32 mm, or 24 x 40 mm) by Kimwipes
- 2. Peel off the protective liner from the iSpacers
- 3. Apply the spacer, with the exposed, tacky side down, onto the dry surface of the slide or coverslip
- 4. Press gently to seal
- 5. To increase the sample well depth, add multiple spacers
- 6. Place the specimen within the well
- 7. Add appropriate amount of RapiClear® into the sample well
- 8. Seal the sample well by placing a coverslip
- Press gently around the edges of the coverslip to ensure a safety seal
- 10. Carefully remove the exceeding solution from the well with Kimwipes
- 11. The area outside the well can be sealed with clear nail-polish to form a hard film

Double-sided sticky iSpacers

- 1. Clean the microscope slide and coverslip (size: 22 x 22 mm, 24 x 32 mm, 24 x 40 mm, or 24 x 60 mm) by Kimwipes
- 2. Peel off the protective liner from the iSpacers
- 3. Apply the spacer, with the exposed, tacky side down, onto the dry surface of the slide or coverslip
- 4. Press gently to seal
- 5. To increase the sample well depth, add multiple spacers
- 6. Place the specimen within the well
- 7. Add appropriate amount of RapiClear® into the sample well
- Peel off the other liner and seal the sample well by placing a coverslip
- 9. Press gently around the edges of the coverslip to ensure a safety seal
- 10. Carefully remove the exceeding solution from the well with **Kimwipes**

LAYOUT



^{*}one pack of #IS201 or #IS211 inside; ** one pack of #IS204 or #IS214 inside.

Demonstration

STORAGE

Avoid compression and store at room temperature (cannot autoclave). iSpacers are not compatible with organic solvents.

For research use only. Not for use in diagnostic procedures

PRODUCT LIST

Current prices may be obtained from our website (www.sunjinlab.com) or from our distributor.

Cat. #	Number of Well	Well Size and Depth	Volume per Well	Quantity
		(mm)	(µl)	(piece/pack)
Single-sided sticky				
IS101	1	26 x 16, <u>0.15 deep</u>	75	50
IS111 (value pack) IS001	_		400	250 50
IS010 (value pack)	1	26 x 16, <u>0.2 deep</u>	100	250
IS002 IS011 (value pack)	1	26 x 16, <u>0.5 deep</u>	250	50 250
IS003	1	26 x 16, <u>1.0 deep</u>	500	50
IS012 (value pack) IS301		-		250 50
IS311 (value pack)	1	Φ16, <u>0.15 deep</u>	40	250
IS302 IS312 (value pack)	1	Ф16, <u>0.2 deep</u>	50	50 250
IS303	1	Ф16 0 25 J	90	50
IS313 (value pack)	1	Φ16, <u>0.35 deep</u>	90	250
IS304 IS314 (value pack)	1	Ф16, <u>0.55 deep</u>	145	50 250
IS006	4	Φ7, <u>0.15 deep</u>	12	50
IS015 (value pack) IS007				250 50
IS016 (value pack)	4	Φ7, <u>0.2 deep</u>	15	250
IS008 IS017 (value pack)	4	Φ7, <u>0.5 deep</u>	35	50 250
IS009	4	Ф7 1.0 door	75	50
IS018 (value pack)	4	Φ7, <u>1.0 deep</u>	13	250
Double-sided sticky				
IS201 IS211 (value pack)	1	26 x 16, <u>0.05 deep</u>	25	50 250
IS203	1	26 - 16 0 25 days	125	50
IS213 (value pack) IS004	1	26 x 16, <u>0.25 deep</u>	123	250 10
IS013 (value pack)	1	26 x 16, <u>3.0 deep</u>	1,500	100
IS005	1	26 x 16, <u>7.0 deep</u>	3,500	10
IS014 (value pack) IS4010		•		50 50
IS4110 (value pack)	1	54 x 16, <u>0.1 deep</u>	90	250
IS4020 IS4120 (value pack)	1	54 x 16, <u>0.2 deep</u>	180	50 250
IS4025	1	54 x 16, <u>0.25 deep</u>	225	50
IS4125 (value pack) IS4030	1	54 х 10, <u>0.25 цеер</u>	223	250 50
IS4130 (value pack)	1	54 x 16, <u>0.3 deep</u>	270	250
IS4040	1	54 x 16, <u>0.4 deep</u>	360	50
IS4140 (value pack) IS305		-		250 50
IS315 (value pack)	1	Ф16, <u>0.05 deep</u>	13	250
IS306 IS316 (value pack)	1	Φ16, <u>0.1 deep</u>	26	50 250
IS307	1	Ф16, <u>0.3 deep</u>	80	50
IS317 (value pack) IS308		Ψ10, <u>0.3 deep</u>	80	250 50
IS318 (value pack)	1	Φ16, <u>0.25 deep</u>	65	250
IS309	1	Φ16, <u>0.4 deep</u>	110	50
IS319 (value pack) IS330	1	-	700	250 10
IS331 (value pack)		Ф16, <u>3.0 deep</u>	780	100
IS204 IS214 (value pack)	4	Φ7, <u>0.05 deep</u>	4	50 250
IS206	4	Ф7 0.25 daen	18	50
IS216 (value pack)	4	Φ7, <u>0.25 deep</u>	18	250